

02-11-2018

**Amendment No. 1****Sub: Amendment to the referred tender enquiry****Ref.: Tender Enquiry: HITES/PCD/MP/07/CSSD/18-19 Dated: 11.10.2018**

The following changes are being incorporated in the above referred Tender Enquiry Document

**SECTION I****NOTICE INVITING TENDER (NIT)****Existing:****Tender Timeline:**

| Sl. No. | Description  | Schedule          |
|---------|--|-------------------|
| c.      | Closing date & time for submission of online bids  | 05.11.2018, 13:00 |
| d.      | Closing date & time for submission of <b>tender processing fee and EMD in physical form*</b> | 05.11.2018, 14:00 |
| e.      | Time and date of opening of online bids  | 05.11.2018, 14:30 |

**Read As:****Tender Timeline:**

| Sl. No. | Description  | Schedule                 |
|---------|--|--------------------------|
| c.      | Closing date & time for submission of online bids  | <b>21.11.2018, 13:00</b> |
| d.      | Closing date & time for submission of <b>tender processing fee and EMD in physical form*</b> | <b>21.11.2018, 14:00</b> |
| e.      | Time and date of opening of online bids  | <b>21.11.2018, 14:30</b> |

**SECTION-VII****TECHNICAL SPECIFICATIONS**

| <b>Schedule No: 01:-CSSD</b> |  |   |
|------------------------------|--|---|
| <b>Para</b>                  | <b>Existing Specification:</b>   | <b>Read as:-</b>  |
|                              | <b>WASHER DISINFECTOR (300 - 350 L) With Accessories</b>   | <b>WASHER DISINFECTOR (280 - 350 L) With Accessories</b>  |
| page no. 49<br>Para 1 d      | <b>Washer Disinfector</b><br>The washer disinfector shall be supplied with universal rack, <b>6 level racks</b> for instrument tray, full size instrument tray as well as stop valves, anti-suction device and plastic water trap. | The washer disinfector shall be supplied with universal rack for instrument tray, full size instrument tray as well as stop valves, anti-suction device and plastic water trap  |
| page 49 para<br>5            | The washer disinfector shall be supplied with universal rack, <b>6 level racks</b> for instrument tray, full size instrument tray as well as stop valves, anti-suction device and plastic water trap.                              | The washer disinfector shall be supplied with universal rack for instrument tray, full size instrument tray as well as stop valves, anti-suction device and plastic water trap. |
| Page No.50,<br>Para 4.A      | <b>Ultrasonic cleaner</b><br>The units should be a compact free-standing bench model, with a built-in tank manufactured from high-   | The units should be a compact free-standing bench model, with a built-in tank manufactured from high-quality (316/304) stainless steel and a solid-state generator that         |

|                          |  |  |
|--------------------------|--|--|
|                          | quality (316/304) stainless steel and a solid-state generator that sends ultrasonic (approx 40 KHz)  | sends ultrasonic ( <b>approx 38-45KHz</b> )  |
| Page No.51,<br>Para 7.D  | <b>Vacuum Cleaner:</b><br>Should have 60 liters per second air flow, 17 kilopascals suction power.   | Should have <b>50-60</b> liters per second air flow, 17 kilopascals suction power.   |
| Page No.55,<br>Para .b.1 | <b>Horizontal Sterilizer 500 litre</b><br>Construction:<br><b>Chamber &amp; Doors:</b> The chamber should be jacketed to ensure the temperature uniformity in chamber with multiple heaters. Specify the number of heaters.  | <b>Horizontal Sterilizer 500 litre</b><br>Construction: (Chamber & Doors) :<br>The chamber should be jacketed to ensure the temperature uniformity in chamber with multiple heaters. Specify the number of heaters ( <b>in steam generators</b> ).   |
| page no 57<br>para B.1   | <b>Sterilizer 250 L or more and with minimum 4 STU capacity with Accessories:</b><br>Construction: (Chamber & Doors) :<br>The chamber should be jacketed to ensure the temperature uniformity in chamber with multiple heaters. Specify the number of heaters.   | Sterilizer 250 L or more and with minimum 4 STU capacity with Accessories:<br>Construction: (Chamber & Doors) :<br>The chamber should be jacketed to ensure the temperature uniformity in chamber with multiple heaters. Specify the number of heaters ( <b>in steam generators</b> )."                  |
| Page no. 59<br>para 3.7  | <b>Table Top Sterilizer with Accessories:</b><br>Chamber: Should be made of S.S.316 & should comply the Pressure Equipment Directive (PED) &EN 13445 norms. Chamber should have working pressure 2.2 bar & design pressure upto 3.8 bar. Chamber should be equipped with electrically heated jacket for preheating on standby mode | Chamber: Should be made of S.S.316 & should comply the Pressure Equipment Directive (PED) &EN 13445 norms. Chamber should have working pressure <b>2.2 -2.5 bar</b> & design pressure upto <b>3.8-4 bar</b> . Chamber should be equipped with electrically heated jacket for preheating on standby mode. |
| Pag no 60<br>para 11     | <b>Table Top Sterilizer with Accessories:</b><br>11) Water Storage Tank: Sterilizer should have inbuilt water reservoir with storage capacity up to 5 L. The water reservoirs should have easy access for cleaning & to avoid bio film.  | Water Storage Tank: Sterilizer should have inbuilt water reservoir with storage capacity up to <b>3.5-5 L</b> . The water reservoirs should have easy access for cleaning & to avoid bio film.   |
| Page No.64,<br>Para 16.E | <b>Drying Cabinet:</b><br>Approximate Dimension: 600X600X 600 mm   | Approximate Dimension: 600X600X 600 mm or <b>capacity : 220 L or more</b>  |
|                          | <b>OTHER EQUIPMENTS</b>  |  |
| Page 69<br>Para n        | RO should of Eureka Forbes/Ion Exchange / Millipore / Kent / Aquacare / Rions / Water expert make  | RO should of Eureka Forbes/Ion Exchange / Millipore / Kent / Aquacare / Rions / Water expert make/ <b>Delta pure</b>   |
|                          |  | <b>Added Para:- 5% tolerance in dimension for All furnitures</b>   |

**All other contents of the tender enquiry including terms & conditions remain unaltered.**

**Note:**

- i. **Prospective Bidders are also advised to check the website regularly prior to the closing date and time of online submission of bids.**
- ii. **If EMD is submitted in the form of BG, then the validity of the BG should be at least 165 days from the date of tender opening.**